



Parainfluenza-Virus Typ 1/2/3-Antikörper

Testparameter anti-Parainfluenza-Virus 1/2/3-IgA (Serum/Liquor)
anti-Parainfluenza-Virus 1/2/3-IgG (Serum/Liquor)

Material Serum, 1 mL
Liquor, 1 mL

Entscheidungsgrenze	Serum	Negativ	Grauzone
	anti-Parainfluenza-Virus 1/2/3-IgA	< 10 U/mL	10 – 15
anti-Parainfluenza-Virus 1/2/3-IgG	< 10 U/mL	10 – 15	
Liquor			
anti-Parainfluenza-Virus 1/2/3-IgA/IgG	< 10 U/mL		
ASI	< 2		

Methode Elisa

Qualitätskontrolle Zertifikat

Anforderungsschein Download und Analysenposition

Auskünfte Infektionsimmunologie