



## Picorna-Viren-Antikörper

<b>Material</b>	<b>Serum</b> , 1 mL	anti-Picorna-Viren-IgG
		anti-Picorna-Viren-IgA
		anti-Picorna-Viren-IgM
	<b>Liquor</b> , 1 mL	anti-Picorna-Viren-IgG
		anti-Picorna-Viren-IgA
		anti-Picorna-Viren-IgM
		<u>Antikörper-Spezifität-Index</u>

Entscheidungsgrenze	Serum	Negativ	Grauzone
	anti-Picorna-Viren-IgG	< 11 U/mL	11 - 14
	anti-Picorna-Viren-IgA	< 10 U/mL	10 - 15
	anti-Picorna-Viren-IgM	< 10 U/mL	10 - 15
	<b>Liquor</b>		
	anti-Picorna-Viren-IgG	< 11 U/mL	
	anti-Picorna-Viren-IgA	< 10 U/mL	
	anti-Picorna-Viren-IgM	< 10 U/mL	
	ASI	< 2	

<b>Methode</b>	<u>Elisa</u>
<b>Qualitätskontrolle</b>	<u>Zertifikat</u>
<b>Anforderungsschein</b>	<u>Download</u> und <u>Analysenposition</u>
<b>Auskünfte</b>	<u>Infektionsimmunologie</u>